#### miRNAs are promising novel biomarkers

#### Inserm/University Paul Sabtier

NOvember 22, 2011

Niels Montano Frandsen, PhD Product Manager

nmf@exiqon.com







#### Exiqon at a glance

#### Highlights

- Proprietary LNA<sup>™</sup> detection technology
- Unique IP to novel group of miRNA biomarkers
- Established one-stop shop for Life Science products
- Unique platform for Molecular Diagnostics
  products
- Multiple license agreements with world leading companies
- Listed on NASDAQ OMX, Copenhagen ("EXQ")
- 75 employees
- Founded 1996

#### Locations







#### Exiqon at a glance

#### **Business divisions**



Exiqon Life Sciences combines leading-edge scientific expertise in gene expression with our proprietary LNA<sup>™</sup> technology. Our products, services and scientific staff enable life science researchers to make groundbreaking discoveries



Exiqon Diagnostics is the leader in providing technologies for miRNA biomarker detection. Exiqon Diagnostics is dedicated *in collaboration with partners* to develop novel molecular diagnostic tests for early detection of diseases and knowledge based treatment selection.



#### MicroRNA research: Broad product portfolio





### Number of miRNA publications strongly growing











- How to study microRNAs introduction to LNA<sup>™</sup>
- miRCURY LNA<sup>™</sup> microRNA arrays
- Early diagnosis of cutaneous T-cell lyphoma
- Detection of microRNAs in clinical samples and diagnostic development using LNA<sup>™</sup> Universal RT microRNA PCR

........

• Early detection of colorectal cancer from human patient blood plasma





#### Six facts about microRNAs

#### Six Facts

- 1. Short non-coding RNA molecules of 19-22 nucleotides
- 2. Post-transcriptional regulators of mRNA
- 3. 1424 annotated human microRNAs\* (1000 2000 predicted)
- 4. Regulate at least one third of all human genes
- 5. Phylogenetically well conserved
- 6. Altered microRNA expression profiles are associated with a number of diseases (cancer, diabetes, neurological disease, viral infection)

\*Sanger miRBase release 17.0, April 2011 http://www.mirbase.org

#### microRNA Mode of Action



#### Binding to 3' UTR of mRNA

- Imperfect complementarity translational repression mRNA decay
- Perfect complementarity **mRNA cleavage**



#### Reasons to use miRNAs as biomarkers

#### Key roles in pathway regulation & tissue differentiation

- Important regulatory roles in many diseases
- •Regulate expression of key proteins in drug response pathways
- •Predict response to targeted therapies
- Actively secreted into the circulation
- •Relatively small number of genes to profile
- •Huge dynamic range (0 to 40,000 molecules per cell)

#### Highly stable in clinical sample preparations

- archival FFPE material
- serum and plasma (routinely collected, minimally invasive)
- other body fluids

vember 23,



### Analyzing miRNAs – challenges using traditional DNA technology

Feature	Challenge
Very short sequences ~ 22nt The length of a normal capture probe or PCR primer	Hard to achieve sufficient sensitivity - especially with AT-rich sequences. <b>No design possible!</b>
There are many microRNA families, members differ by single nucleotides	Hard to achieve sufficient specificity to discriminate between family members
Highly variable GC content – from 5-95%	A serious problem with multiplex assays (such as microarrays – full length DNA probes have a Tm span of 45 C!)





#### LNA<sup>™</sup> technology overview







#### LNA<sup>™</sup> technology enables design of short RNA analysis tools

The power of LNA™ to modulate affinity of probes for RNA target¹							
DNA oligo	caacatcagtctgataagct	47 C					
LNA/DNA mixmer (1 LNA)	caacatcagTctgataagct	52 C					
LNA/DNA mixmer (2 LNAs)	caacatcagTcTgataagct	62 C					
LNA/DNA mixmer (4 LNAs)	caacaTcagTcTgaTaagct	73 C					
LNA/DNA mixmer (6 LNAs)	caaCaTcagTcTgaTaAgct	78 C					
LNA/DNA mixmer (8 LNAs)	cAaCaTcagTcTgaTaAgCt	84 C					

Discriminative power of LNA<sup>™ 2</sup>

Target Probe	Perfect match DNA	Single mismatch 3'ccaggaa <mark>g</mark> gaaccac-5'	Δ <b>T</b> <sub>m</sub>
DNA 15mer 5'ggtccttacttggtg3'	T <sub>m</sub> = 59,4 C	T <sub>m</sub> =55,7 C	3,7 C
DNA/LNA 15mer 5'ggtccttActtggtg3'	T <sub>m</sub> =60,9 C	T <sub>m</sub> = 52,7 C	8,2 C

- LNA enhances <u>affinity</u> of oligonucleotides for complementary RNA and DNA targets
- LNA<sup>™</sup> technology enables <u>probe design</u> with a desired Tm for short RNA targets

 LNA enhances <u>specificity</u> of oligonucleotides for complementary RNA and DNA targets

<sup>1</sup> Capital letters in red indicate LNA<sup>™</sup>

<sup>2</sup>  $\Delta$ Tm is the difference between Tm of perfect match and mis-match duplexes. Capital letters indicate LNA<sup>TM</sup>

11



#### Sensitivity and specificity of LNA<sup>™</sup> (Locked Nucleic Acid) technology



In situ hybridizations using DIG-labeled probes against liver specific and highly expressed miR-122.

ovember 23, 2011

Specificity achieved with LNA



s - somites, h - heart

Pictures kindly provided by D. Sweetman, Münsterberg Group, University of East Anglia, Norwich, UK

miR-206 UGGAAUGUAAGGAAGUGUGUGG miR-1 UGGAAUGUAAAGAAGUAUGUAU

12

. .



### T<sub>m</sub> normalization of microRNA inhibitors achieved by LNA design

# Tm normalization ensures high potency of LNA inhibitors regardless of target microRNA GC content



- How to study microRNAs introduction to LNA<sup>™</sup>
- miRCURY LNA™ microRNA arrays
- Early diagnosis of cutaneous T-cell lyphoma
- Detection of microRNAs in clinical samples and diagnostic development using LNA<sup>™</sup> Universal RT microRNA PCR

........

• Early detection of colorectal cancer from human patient blood plasma





#### DNA capture probes have poor sensitivity for low Tm targets and decreased specificity for high Tm targets





### Development of miRCURY LNA<sup>™</sup> microRNA array platform

Capture probes are experimentally validated





#### Testing specificity of capture probes with synthetic miRNA pools





#### miRCURY LNA<sup>™</sup> microRNA arrays have Tm normalized capture probes



Using LNA<sup>™</sup>, the Tm is increased significantly and the Tm range is narrowed significantly, compared to DNA probes. This results in increased stringency and optimal hybridization conditions for the LNA<sup>™</sup> capture probes.



#### Sensitivity of 6<sup>th</sup> generation miRCURY LNA<sup>™</sup> microRNA arrays





## miRCURY LNA<sup>™</sup> microRNA Arrays have unmatched sensitivity for ALL microRNAs – even the AT-rich



LNA™ capture probes

DNA capture probes



% GC-content in the microRNAs

DNA-based arrays fail to detect low GC-content microRNAs (>20% of all microRNAs)





## Excellent discrimination between closely related microRNAs – let-7 family specificity

Array capture probes

#### hsa-let-7 family sequences

hsa	let-7a	UGA	GGU	AGU	AGG	UUG	UAU	AGU	U
hsa	let-7b	UGA	GGU	AGU	AGG	UUG	U <mark>G</mark> U	<mark>G</mark> GU	U
hsa	let-7c	UGA	GGU	AGU	AGG	UUG	UAU	<mark>G</mark> GU	U
hsa	let-7d	<mark>A</mark> GA	GGU	AGU	AGG	UUG	<mark>C</mark> AU	AGU	-
hsa	let-7e	UGA	GGU	AG <mark>G</mark>	AGG	UUG	UAU	AGU	-
hsa	let-7f	UGA	GGU	AGU	AG <mark>A</mark>	UUG	UAU	AGU	U
hsa	let-7g	UGA	GGU	AGU	AG <mark>U</mark>	UUG	UA <mark>C</mark>	AGU	-
hsa	let-7i	UGA	GGU	AGU	AG <mark>U</mark>	UUG	U <mark>GC</mark>	<mark>U</mark> GU	-

#### microRNA spike-in

	let-7a	let-7b	let-7c	let-7d	let-7e	let-7f	let-7g	let-7i
let-7a	100	2	17	4	4	2	1	2
let-7b	1	100	4	1	1	1	1	1
let-7c	0	8	100	0	1	0	0	0
let-7d	2	2	5	100	1	0	0	0
let-7e	1	0	0	0	100	0	0	0
let-7f	6	3	5	3	3	100	2	3
let-7g	0	0	1	0	0	1	100	4
let-7i	0	3	0	0	0	0	2	100

21



#### Coverage – miRBase v. 16.0





#### FFPE samples –

### testing the performance of the miRCURY LNA™ microRNA Array





#### microRNAs from FFPE vs. fresh frozen tissue





## Excellent correlation of microRNA profiles from FFPE and fresh frozen samples





## Minimal batch to batch variation – ideal for single color expression profiling



Very low CV between different batches of the array:

median CV  $\leq 10\%$ 

median CV 2-4% typically obtained between different batches

enabling superior single color experiments



- How to study microRNAs introduction to LNA<sup>™</sup>
- miRCURY LNA<sup>™</sup> microRNA arrays
- Early diagnosis of cutaneous T-cell lyphoma
- Detection of microRNAs in clinical samples and diagnostic development using LNA<sup>™</sup> Universal RT microRNA PCR

........

• Early detection of colorectal cancer from human patient blood plasma





- CTCL are the most frequent primary lymphomas of the skin
- Early diagnosis of CTCL has important consequences for therapeutic options and determination of prognosis
- Early diagnosis of CTCL is difficult due to the great clinical, pathological, and histological resemblance to benign inflammatory skin diseases (psoriasis and eczema).
- From initial skin lesion to definite CTCL diagnosis: 6 years (median)



\* Farber EM, Nall L. Epidemiology: Natural History and Genetics. In: Roenigk HH Jr, et al, editors. Psoriasis, 3<sup>rd</sup> edn. Marcel Dekker, New York, 1998:107-157.



#### Sample selection and experimental set up

		Cutane	eous Lympho	oma ( <i>n</i> =63)		Benign skin disease or normal skin ( <i>n</i> =85)			Durahua			
	MF	SS	CALCL	NOS	Σ	AD	ND	PP	PN	NN	Σ	P-value
Age (years) *												<0.001
<30	0	0	1	0	1	19	0	4	4	2	29	
30-44	5	0	0	0	5	1	0	6	3	0	10	
45-59	9	1	1	2	13	0	1	19	7	0	27	
60-74	14	6	2	2	24	0	3	12	2	0	17	
≥75	10	0	2	5	17	0	0	1	1	0	2	
Gender *												1.00
Male	23	7	5	5	40	8	1	30	16	2	57	
Female	15	0	1	4	20	12	3	12	1	0	28	
Microarray batch												1.00
1	8	2	1	2	13	4	1	8	3	1	17	
2	7	2	2	1	12	4	1	8	3	1	17	
3	7	1	2	2	12	4	1	8	3	0	16	
4	9	1	2	2	14	4	0	9	4	0	17	
5	8	1	1	2	12	4	1	9	4	0	18	

• For each sample RNA was extracted from 6 x 10µM FFPE sections (archival).

• Samples were divided into a training set of 90 and a test set of 58 (randomized) samples.

• 100ng total RNA was used per array.

Reference: Ralfkiaer et al, 2011. Blood doi:10.1182/blood-2011-06-358382



#### EXIQON LNA™ microRNA arrays identify 27 highly significant microRNAsedifferentially expressed between CTCL and benign skin disease



Reference: Ralfkiaer et al, 2011. Blood doi:10.1182/blood-2011-06-358382

30

# Identification of a microRNA classifier for CTCL versus benign skin disease



The 5 most differentially expressed microRNAs from the microarray data make up the classifier: miR-203, miR-205, miR-326, miR-663b and miR-711.

The classifier is validated in the test sample set and in two independent sample sets including a mouse xenograph model.



EXIQON Seek Find Verify

doi:10.1182/blood-2011-06-358382

31

#### Validation of a highly sensitive and specific microRNA classifier



miR-155 functions as an oncogene in other hematological malignancies

miR-203 is important in keratinocyte development and may act as a tumor suppressor

miR-205 is a tumor suppressor possibly through targeting VEGF



From bloodjournal.hematologylibrary.org at HOUSTON ACADEMY OF MEDICINE on September 7, 2011. For personal use only.



Prepublished online August 24, 2011; doi:10.1182/blood-2011-06-358382

#### Diagnostic microRNA profiling in cutaneous T-cell lymphoma (CTCL)

Ulrik Ralfkiaer, Peter H. Hagedorn, Nannie Bangsgaard, Marianne B. Løvendorf, Charlotte B. Ahler, Lars Svensson, Katharina L. Kopp, Marie T. Vennegaard, Britt Lauenborg, John R. Zibert, Thorbjørn Krejsgaard, Charlotte M. Bonefeld, Rolf Søkilde, Lise M. Gjerdrum, Tord Labuda, Anne-Merete Mathiesen, Kirsten Grønbaek, Mariusz A. Wasik, Malgorzata Sokolowska-Wojdylo, Catherine Queille-Roussel, Robert Gniadecki, Elisabeth Ralfkiaer, Carsten Geisler, Thomas Litman, Anders Woetmann, Christian Glue, Mads A. Røpke, Lone Skov and Niels Odum

- How to study microRNAs introduction to LNA<sup>™</sup>
- miRCURY LNA<sup>™</sup> microRNA arrays
- Early diagnosis of cutaneous T-cell lyphoma
- Detection of microRNAs in clinical samples and diagnostic development using LNA™ Universal RT microRNA PCR

........

• Early detection of colorectal cancer from human patient blood plasma





### Circulating microRNAs in cell-free serum and plasma

- Major potential as biomarkers

lung, colon cancer and diabetes have specific serum-miRNA profiles	Chen <i>et al.,</i> 2008, Cell Res, 1-10
MiR-92 is significantly elevated in plasma of patients with CRC, potential non-invasive molecular marker for CRC screening	Ng <i>et al.</i> , 2009, Gut, <b>58:</b> 1375- 1381
Serum miR-21 associated with relapse-free survival in diffuse large cell B <b>lymphoma</b>	Lawrie <i>et al.,</i> BJH, 2008
human prostate cell line xenotransplanted into mice and tumor-derived human specific miRNAs identified in plasma	Mitchell <i>et al.,</i> 2008, PNAS, 105(30):10513.
microRNAs differentially expressed in serum of ovarian cancer patients	Resnick <i>et al.,</i> 2009 Gyn Oncology 112: 55-59
microRNAs differentially expressed in circulating tumor-derived <b>exosomes</b> of <b>ovarian cancer</b> patients	Taylor <i>et al.,</i> 2008, Gyn Oncology, 110:13
The circulating microRNA-221 level in patients with malignant melanoma as a new tumor marker	Kanemaru H <i>et al</i> ., 2011, J Dermatol Sci. 61(3):187-93
Plasma miR-208 as a Biomarker of Myocardial Injury	Ji <i>et al.,</i> 2009, Clin. Chem. 55: 11
Plasma miRNA signature for diabetes mellitus	Zampetaki <i>et al.,</i> 2010, Circ Res. 2010:107:810-817



#### Clinical QC of miRNA stability in blood derived Plasma

#### UNIVERSITY OF COPENHAGEN

#### FACULTY OF LIFE SCIENCES

### Stability of miRNAs in plasma stored at room temperature before extraction.



Place, date, unit, occasion etc. Slide 1

36


# Clinical QC of miRNA stability in blood derived Plasma





# microRNAs as a source of biomarkers in blood



Figure 1 | Sources of blood-based biomarkers. Interactions between cancer cells, tumor microenvironment and the host result in the release of circulating cellular and molecular elements that can serve as rich sources of circulating biomarkers.

S Hanash et al Nat Rev Clin Oncol. 2011 Mar;8(3):142-50.





39

# microRNAs as a source of biomarkers in blood



## REVIEWS MicroRNAs in body fluids—the mix of hormones and biomarkers Maria Angelica Cortez, Carlos Bueso-Ramos, Jana Ferdin, Gabriel Lopez-Berestein, Antl K. Sood and George A. Calin

7248–7259 Nucleic Acids Research, 2010, Vol. 38, No. 20 doi:10.1093/nar/gkq601 Published online 7 July 2010

# Export of microRNAs and microRNA-protective protein by mammalian cells

Kai Wang, Shile Zhang, Jessica Weber, David Baxter and David J. Galas\* $^{\star,\dagger}$ 

# Argonaute2 complexes carry a population of circulating microRNAs independent of vesicles in human plasma

Jason D. Arroyo<sup>®</sup>, John R. Chevillet<sup>®</sup>, Evan M. Kroh<sup>®</sup>, Ingrid K. Ruf<sup>®</sup>, Colin C. Pritchard<sup>®</sup>, Donald F. Gibson<sup>®</sup>, Patrick S. Mitchell<sup>®,1</sup>, Christopher F. Bennett<sup>&,c</sup>, Era L. Pogosova-Agadjanyan<sup>d</sup>, Derek L. Stirewalt<sup>d</sup>, Jonathan F. Tait<sup>®</sup>, and Muneesh Tewari<sup>A,&a</sup>z

cell biology

MicroRNAs are transported in plasma and delivered to recipient cells by high-density lipoproteins Kasey C. Vickers<sup>12</sup>, Brian T. Palmisano<sup>1</sup>, Bassem M. Shoucri<sup>1</sup>, Robert D. Shamburek<sup>1</sup> and Alan T. Remaley<sup>1</sup>

Cortez, M. A. et al. Nat. Rev. Clin. Oncol. 8, 467-477 (2011)





# Limited clinical samples





# miRCURY LNA<sup>™</sup> Universal RT microRNA PCR - product overview

## Ready-to-use panels (384-well), for full screening:

- Human, 742 microRNAs, 2 plates: 40 ng total RNA or 70 µl plasma/serum
- Mouse/Rat, 751 microRNAs, 2 plates: 40 ng total RNA or 70 μl plasma/serum

## Custom Ready-to-use panels (96- and 384-well):

Design your own Pick&Mix Panel – include only your microRNAs of interest

## Individual primer sets, pre-designed and validated (human, mouse, rat)

- 1250 human, mouse and rat microRNA assays
- Reference gene primer sets (small RNAs)
- 1 pg total RNA is sufficient for reliable quantification

## Custom LNA<sup>™</sup> primers sets – design on-line (complex algorithm)

## Reagents – optimized for use with microRNA primer sets:

- cDNA synthesis kit (86-32 rxns)
- SYBR<sup>®</sup> Green master mix (2.5 ml and 25 ml)

## Exiqon's GenEx data analysis software















# The challenge of miRNA qPCR:



# Strategies for microRNA qPCR: PolyA tailing, Universal RT, SYBR® Green detection



## Advantages:

- One RT reaction
- Melt curve analysis is possible

## **Disadvantages:**

- Poor specificity: DNA primer close to length of microRNA = high background and poor sensitivity
- Many AT rich sequences will not be detectable
- Mismatches in the middle and 5' end difficult to discriminate

# Best of two worlds: Universal RT with two microRNA specific PCR primers

Step 1: First-strand synthesis (RT)



# Polyadenylation

Reverse transcription

# Step 2: Real-time PCR amplification



Two LNA<sup>™</sup> enhanced microRNA specific primers

SYBR Green detection

# Simple and rapid system for genome profiling of miRNA



## **1 RT reaction per sample**

- Minimal sample input
- Few handling steps

## 1 one reaction per well

- Small reaction volume 10 μl
- No assay pipetting just add cDNA and SYBR® Green master mix
- **3 hours protocol**

No need for pre-amplification

Plate layout file for direct import to GenEx



# miRCURY LNA™ Universal RT microRNA PCR System





# LNA<sup>™</sup>-enhanced primers enable detection of ALL microRNAs – even the AT-rich





# Excellent family member specificity

hsa-let-7a	UGAGGUAGUAGGUUGUAUAGUU
hsa-let-7b	UGAGGUAGUAGGUUGU <mark>G</mark> UGGUU
hsa-let-7c	UGAGGUAGUAGGUUGUAU <mark>G</mark> GUU
hsa-let-7d	AGAGGUAGUAGGUUG <mark>C</mark> AUAGU-
hsa-let-7e	UGAGGUAG <mark>G</mark> AGGUUGUAUAGU-
hsa-let-7f	UGAGGUAGUAG <mark>A</mark> UUGUAUAGUU
hsa-let-7g	UGAGGUAGUAG <mark>U</mark> UUGUA <mark>G</mark> AGU-
hsa-let-7i	ugagguaguag <mark>u</mark> uugu <mark>ggu</mark> gu-

		hsa-let-7a	hsa-let-7b	hsa-let-7c	hsa-let-7d	hsa-let-7e	hsa-let-7f	hsa-let-7g	hsa-let-7i
ets	hsa-let-7a	100%	0.0%	5.2%	0.0%	0.0%	0.1%	5.1%	0.0%
r S(	hsa-let-7b	0.0%	100%	0.5%	0.0%	0.0%	0.0%	0.0%	0.0%
ne	hsa-let-7c	1.2%	0.0%	100%	0.0%	0.0%	0.0%	0.1%	0.0%
orir	hsa-let-7d	0.0%	0.0%	0.0%	100%	0.0%	0.0%	0.0%	0.0%
Ā	hsa-let-7e	0.4%	0.0%	0.1%	0.3%	100%	0.0%	0.0%	0.0%
RN	hsa-let-7f	0.0%	0.0%	0.0%	0.0%	0.0%	100%	0.0%	0.0%
cro	hsa-let-7g	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	100%	0.0%
ä	hsa-let-7i	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	100%

# Let-7 family templates



# Blood collection tubes and subsequent effect on microRNA qPCR







# Correlation between microRNA input and median Ct – near perfect linearity at very low concentrations as found in plasma and serum





# Genome-wide microRNA profiling based on serum

Profiling of microRNAs from 35µl of blood derived serum – with no pre-amplification



# miRNA levels varies between healthy individuals, but within discrete ranges





# **QC Control: 490,368 PCR data point sample integrity check**



## Exiqon Hsa-miR-Databank Profile Master Average Human Serum

complete miRNA serum profile (numerically ordered)

# Serum/plasma focus panel - 175 assays



Profile all relevant serum/plasma microRNA from 20 µl serum/plasma

96-well (1 panel in 2 plates)

384-well(1 plate with 2 panels)

All the relevant controls



# Human Urine microRNA Profile



## Exiqon Hsa-miR Databank Profile. Master Average Human Urine

Complete miRNA Urine Profile (Numerically ordered)



- How to study microRNAs introduction to LNA<sup>™</sup>
- miRCURY LNA<sup>™</sup> microRNA arrays
- Early diagnosis of cutaneous T-cell lyphoma
- Detection of microRNAs in clinical samples and diagnostic development using LNA<sup>™</sup> Universal RT microRNA PCR

........

Early detection of colorectal cancer from human patient blood plasma





# **Detection of early-to-late stage pancreas cancer in blood plasma**

## The experiment

- 30 pancreas cancer:
  - 3 local
  - 19 locally advanced
  - 8 metastatic
- 6 chronic pancreatitis
- 14 normal pancreas
- Plasma samples
- Whole genome microRNA profiling using LNA Universal RT microRNA PCR Panels
- Unsupervised clustering based on all detected miRNAs
- Normalized data
- Pancreatitis plasma clusters with normal plasma, and away from cancer plasma, irrespective of cancer stage





# **Detection of early-to-late stage pancreas cancer in plasma**

## Cancer vs. normal in plasma

- Assumptions:
  - Pancreatitis and normal grouped as normal
  - All cancer stages grouped as cancer
- Clustering based on most significant miRNAs
- Promising separation between malignant and nonmalignant pancreatic disease in plasma





## **Case study**



# microRNA profiles can discriminate tumor type in serum samples





# Unmet need for early detection test of Colorectal Cancer

## **Colorectal cancer**

- Estimated new cases in 2010 in the USA:
  - Colon: 102,900
  - Rectum: 39,670
- Estimated deaths in 2010 in the USA:
  - Combined: 51,370

## • Early diagnosis results in resectable cancer with much improved prognosis





# Rationale for early detection test of CRC in blood plasma

## **Colorectal cancer (CRC) screening guidelines**

For individuals between 50 and 75:

Colonoscopy every 5 years (False-negative rate ~5%)

## or

- Annual FOBT (False-negative rate 20-75%)
- Poor compliance: <50% of population are screened</li>

Large unmet need for minimally invasive screening assay for detection of CRC

# Proposed application of early detection blood test in Colorectal cancer (CRC) screening Symptoms Colonoscopy Age: > 50 years miRNA blood test If positive, then colonoscopy If negative, come back next year





# Work flow compatible with standard clinical procedures

	Total lab time = 5 hours			
	time = 1 hour	time = 1.5 hours	time = 2 hours	time = 0.5 hours
Plasma sample	RNA isolation	RT reaction	PCR reaction	Data QC and analysis.
<ul> <li>Standard Hospital procedures</li> </ul>	<ul> <li>RNA isolation 200ul plasma required</li> <li>Isolation of total RNA for from sample (no bias)</li> </ul>	<ul> <li>Standard RT reaction.</li> <li>Robust and reproducible</li> </ul>	<ul> <li>PCR with LNA specific primers for detection of miRNAs</li> </ul>	<ul> <li>Quality control</li> <li>Data analysis and results</li> <li>Diagnostic interpretation</li> </ul>



- Assessible to standard hospital protocols
- No special handling/ storage requirements
- Low volume requirement



- FDA approved technology
- Conventional qPCR
- No black box algorithms
- Automatable





# Normalization of qPCR data is not straightforward

## Normalization of qPCR data is necessary to make expression values comparable across samples

Normalization adjusts for technical biases:

- variable amount of sample input RNA
- variable sample input RNA quality
- variable assay efficiencies due to e.g. sample specific assay inhibition

## Prior assumptions about housekeeping genes across the samples of interest is not needed





64

# Blank purification as a QC control of miRNA detection

Multiple controls on miRNA panel together with blank púrification produce good QC process

- Every system has a background and this needs to be assessed critically
- Water sample is exposed to the entire process to assess true background and thresholds.
- This is critical for the analysis of miRNAs in serum and other biofluids where levels are very low.





# Extensive QC and data analysis pipeline

## QC pipeline

- Automated, user-guided QC pipeline
- Data flagging and cleanup customized for qPCR platform
- Diverse normalization protocols implemented
- Comprehensive QC report
- Visualization of potential technical and sample biases
- Identify samples that may be affected by contamination from blood cells





66

# **Pre-screen: Profile differs in plasma from CRC patients and healthy controls**

- 50 stage II colorectal cancer patients and 50 age- and sex-matched colonoscopy negative controls
- Plasma samples (pre-endoscopy)
- Screening defined 378 candidate miRNAs present in plasma







# Focused panel of microRNAs in plasma may be used as biomarker for CRC



	All Samples	Hospital 1	Hospitals 2-5
Sensitivity *	75%	80%	73%
Specificity *	80%	78%	82%
(n) Cancer	151	49	102
(n) Control	76	36	40

\* The same cutoff score was applied on all samples in the study



# **Development of miRNA Early Detection Test of CRC in blood plasma**

DISCOVERY PHASE				VALIDATION PHASE
Genome wide screening	Normalization, QC, processing	Candidate miRNA discovery screen.	Bioinformatics, data analysis,	Validation Set miRNA signature .
<ul> <li>50 controls</li> <li>50 CRC patients</li> <li>742 miRNA screen</li> <li>Genome wide</li> </ul>	<ul> <li>Multiple QC check</li> <li>Data flagging</li> <li>Normalization</li> </ul>	<ul> <li>76 controls</li> <li>151 CRC patients</li> <li>Focused Blood panel</li> <li>Multiple controls</li> </ul>	<ul> <li>Data analysis</li> <li>Quality control</li> <li>ROC curve</li> <li>miRNA selection</li> </ul>	<ul> <li>3000 patients (2011)</li> <li>Defined miRNA signature</li> <li>Pick &amp; Mix panel</li> <li>Multiple controls</li> </ul>
			Hospital 1 (n=85)	



69

# miRNA Early Detection Biomarker from Discovery to Development



# **Biomarker discovery workflow and panel selection options**



Serum/plasma Focus microRNA PCR Panel

- No need to spend time and money on miRNome-wide screening (Exiqon has done this for you)
- Comfortably start your focused screening with thoroughly selected serum/plasma microRNAs
- Move to Exigon's Pick&Mix microRNA PCR Panel once microRNAs of interest are identified

# Conclusions and future prospects

- Circulating microRNAs are promising biomarkers
- MicroRNA biomarker discovery should be included in early clinical studies
- microRNAs can be sensitively and reproducibly detected in serum / plasma without pre-amplification using the miRCURY LNA™ Universal RT microRNA PCR system
- Simple workflow allows test results from 100 µl of plasma within one day
- Flexible platform enables whole genome profiling and focused / customizable panels
- A plasma microRNA signature for colorectal cancer was developed and performs well in samples from independent hospitals
- Validation of the signature in a larger cohort is on-going



## **Acknowledgements**

### Exiqon, Denmark & US

Adam Baker Maria Wrang Teilum Christian Glue Stine Jørgensen Thorarinn Blondal Søren Jensby Nielsen Randi Holm Anita Pacht Lene Bach Peter Mouritzen Jannie Christensen Mette Carlsen Mohr Søs Ludvigsen Ditte Andreasen Hanni Willenbrock

## Exiqon Marketing Ina Dahlsveen Roman Kurek Niels Frandsen

**Copenhagen University** Nils Brünner Jan Stenvang

**Hvidovre Hospital** Hans Jørgen Nielsen Ib Jarle Christensen

## **Aarhus University Hospital Skejby**

Torben F. Ørntoft Claus L. Andersen Steffen Grann Jensen

## Funding:

**The Danish Advanced Technology Foundation** Platform grant 007-2009-2

Concerning the miRCURY LNA™ Universal RT microRNA PCR system:

NOTICE TO PURCHASER: LIMITED LICENSE Purchase of this product includes an immunity from suit under patents specified in the product insert to use only the amount purchased for the purchaser's own internal research. No other patent rights are conveyed expressly, by implication, or by estoppel. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA. SYBR Green is a trademark of Invitrogen.
## Thank you for you attention!

miRNAs are promising novel biomarkers

## Inserm/University Paul Sabtier

NOvember 22, 2011

Niels Montano Frandsen, PhD Product Manager

nmf@exiqon.com



